

**Adamantane (<sup>1</sup>H/<sup>13</sup>C Pulse Calibration/Shimming/Chemical Shift Reference)**


**I. <sup>1</sup>H pulse calibration** 2

**II. <sup>13</sup>C pulse calibration** 4

## I. Adamantane – $^1\text{H}$ pulse calibration


**NOTE:** Adamantane has a high-degree of symmetry resulting in a high degree of motion. Thus both  $^{13}\text{C}$  and  $^1\text{H}$  lines will be much narrower than observed in typical organic solids.

- (1) Spin adamantane at 10 kHz and check tuning. Re-tune if necessary
- (2) Using one-pulse experiment, acquire  $^1\text{H}$  spectrum and set adamantane  $^1\text{H}$  to be on-resonance.




**Figure 1:**  $^1\text{H}$  spectrum of adamantane.

**NOTE:** For  $^1\text{H}$  and  $^{13}\text{C}$  pulse calibrations it is important to set recycle delay to at least  $5 \times T_1$ . Insufficiently long delays will shift the maxima of nutation curves (**Figure 2**).



**Figure 2:** Effect of recycle delays with respect to  $T_1$  on nutation curve. From Glenn Facey NMR Blog: <http://u-of-o-nmr-facility.blogspot.com/2009/06/90-degree-pulse-determinations.html>

- (3) To calibrate pulse width – measure nutation curves as a function of  $^1\text{H}$  pulse power with the major adamantane peak on-resonance (**Figure 3**).



**Figure 3:** Nutation curve for  $^1\text{H}$  pulse length of adamantane.


- (4) Calculate the  $90^\circ$  pulse length from the signal null occurring at the  $180^\circ$  or  $360^\circ$  pulse lengths.

**NOTE:** For a typical solids experiment determine the power level for  $2.5\mu\text{s}$  (100kHz) and  $4\mu\text{s}$  (62.5kHz)  $90^\circ$  pulses. The actual power that the probe can handle and amplifiers can provide will be specific to the current setup.

## II. Adamantane – $^{13}\text{C}$ pulse calibration

**NOTE:** Adamantane has a high-degree of symmetry with a high degree of motion resulting in much smaller  $^{13}\text{C}$ - $^1\text{H}$  dipolar coupling than observed in typical organic solids. With moderately fast MAS rates (10 kHz), low power  $^1\text{H}$  decoupling is sufficient to decouple the  $^1\text{H}$ 's. Typically this power is less than 0.1W, but will be dependent on the probe/amplifier setup.


- (1) Spin adamantane at 10 kHz and use a single-pulse with  $^1\text{H}$  decoupling experiment to acquire a  $^{13}\text{C}$  spectrum of adamantane.
- (2) Set the  $^1\text{H}$  decoupling frequency to that determined previously for the  $^1\text{H}$  spectrum. The  $^1\text{H}$  decoupling power will be low ( $\sim 0.1\text{W}$ ) which will allow for long acquisition times (200 to 300ms).
- (3) Set the carrier to be between the two  $^{13}\text{C}$  resonances (**Figure 4**)



**Figure 4:**  $^{13}\text{C}$  spectrum of adamantane, with proper shimming (FWHM  $\sim 4.5$  Hz), collected using ultra low power  $^1\text{H}$  decoupling.

- (4) To calibrate pulse width, measure  $^{13}\text{C}$  nutation curve (**Figure 5**) – peak height as a function of  $^{13}\text{C}$  pulse power. Determine power levels for  $90^\circ$  pulse lengths of  $5\mu\text{s}$  (50kHz), and  $4\mu\text{s}$  (62.5kHz) (depending on power handling of the probe).

**NOTE:** For  $^1\text{H}$  and  $^{13}\text{C}$  pulse calibrations it is important to set recycle delay to at least  $5 \times T_1$ . Insufficiently long delays will shift the maxima of nutation curves (**Figure 2**).




**Figure 5:** Nutation curve for  $^{13}\text{C}$  pulse length on adamantane, collected with  $^{13}\text{C}$  one-pulse with  $^1\text{H}$  decoupling pulse program.

#### (5) Adamantane – Shimming

- At the magic-angle, shims will behave differently than solution samples see (A. Sodickson, D.G. Cory, Shimming a High-Resolution MAS Probe, J. Magn. Reson. (1997) 128: 87-91 <https://doi.org/10.1006/jmre.1997.1218> )
- Align the probe with the spinning module rotor axis along either the X or Y – axis of the RT shims.

**NOTE:** If the rotor-axis is aligned along the Y-axis, applying current to the X shim will minimally affect the  $^{13}\text{C}$  line width.

- One should be able to obtain a FWHM below 10 Hz for the  $^{13}\text{C}$  adamantane resonances. Typically, Z and Y shims (or X depending on direction that probe is oriented) can reduce the  $^{13}\text{C}$  line width. Large values of X2-Y2 will typically reduce the asymmetric foot at the base of the  $^{13}\text{C}$  resonance.



- (6) Adamantane can be used as an external chemical shift standard by setting the downfield (left-most)  $^{13}\text{C}$  resonance to 38.48ppm. (see Morcombe CR, Zilm KW. Chemical shift referencing in MAS solid state NMR. J Magn Reson. 2003 162:479-86. [https://doi.org/10.1016/s1090-7807\(03\)00082-x](https://doi.org/10.1016/s1090-7807(03)00082-x) )